



Why is Topi-Click® a better way to dispense medication than using a tube or jar?

- Topi-Click is very easy to use for dispensing measured doses.
- A patient does not need to use her or his hands to apply the medication, decreasing the risk of their loved ones and pets contacting the drug.
- The full dose is applied to the intended area instead of being left on the hands and wasted when washed off.

TOPI-CLICK® makes it easy to accurately measure the dose and apply topical medications to the intended area, thereby decreasing the risk of improper dosing and contamination while improving compliance and outcomes of therapy.



Watch the video at www.Topi-Click.com to see the Topi-Click in action and the highlights from Oprah, Dr. Phil, and The Doctors national talk shows.

Healthcare Provider Benefits of **TOPI-CLICK®**

- Improves compliance for better outcomes of therapy through ease-of-use and convenience.
- Audible, visual and tactile clicks make it easy for even your most challenged patients to accurately dose medications.
- Demonstration on Oprah, Dr. Phil and The Doctors national talk shows made Topi-Click nationally recognized as the dosing applicator used for topicals.
- Distinguishable and elegant design that's unique from any other type of dispensing system to differentiate your treatments
- Domed applicator pad allows medication to be applied without touching medication with the hands, thereby reducing the risk of spreading medicine to loved ones and pets.
- Easier to use than oral dosing syringes (doesn't resemble a hypodermic needled syringe)
- More pharmaceutically elegant than tubes or jars
- Available in multiple colors (Blue, Pink and more) for differentiating various medications
- Promotes additional compliance with a "Refill" reminder system so there are no gaps in treatments
- User friendly enough to be operated in the dark and almost impossible to discharge accidentally
- Dispenses an average metered dose of 0.25 ml per click.

TOPI-CLICK®

SPECIFICATIONS

Each click (¼ turn of the Topi-Click® container base) moves the plunger (topical elevator) approximately 0.25ml of volume space, thereby dispensing 0.25ml of cream on average. The dispensing mechanism is similar to a syringe, but with more control of the plunger movement, so Topi-Click® is easier to use.

The Topi-Click® container can be used to dispense a total of approximately 5ml to 35ml of a topical preparation.

Topi-Click® is made of FDA medical-grade polypropylene plastic. UV block plastic protects the contents from UV rays.



Topi-Click® can be labeled like a prescription vial.

The domed applicator pad allows patient to apply the medication without dispensing onto hands first.

Audible, visual and tactile clicks make it easy for even your most challenged patients to accurately dose medications.

A “Refill” reminder system helps patients to reorder in a more timely manner.

Maximize the Benefits of Topical Medications

**DOSE IT EASY
AND QUICK**



Topi-Click® Benefits:

- Helps to improve compliance with topical treatments
- Ease-of-use increases patient satisfaction from treatments.
- Has a “Refill” indicator for prompting patients to act so there are no gaps in treatment.
- Demonstrations on Oprah, Dr. Phil and The Doctors national talk shows have made Topi-Click nationally recognized as THE dosing applicator for topical medications.

