



Visit www.Topi-Click.com to view full pages:

[Topi-Click selected as Technology Innovation Award nominee at the Drug Delivery Partnerships conference 2010](#)

In its 14th year, Drug Delivery Partnerships has claimed its place as the leading worldwide event for pharmaceutical, biotech and drug delivery professionals engaged in leading pipeline expansion through drug delivery innovation.

The Topi-Click TDA protects the public and the drugs' success with features that help promote and attain compliance.

Friday, January 22, 2010
2010 DDP Awards: Technology Innovation Nominee: Topi-Click

Technology: The patented Topi-Click topical dosing applicator (TDA)

"Dose it easy and quick with Topi-Click."

"We don't put deodorant on with our hands so why do we ask patients to put topical meds on with their hands?" The patented Topi-Click topical dosing applicator (TDA) provides a pill bottle's ease of use with the dosing precision of syringes to solve this problem. Topi-Click provides an integrated applicator and with a metered dosing clicker to overcome the non-compliance procedure of having to rub topical meds onto treatment areas with the hands. We know the hand is a great absorber of topical; therefore, the patient is using the treatment area to rub a large percentage of the medicine into their hands. Using the hands can also increase chances of others accidentally being exposed to someone else's prescription when holding or shaking hands. Very recently, a male topical testosterone drug was delayed approval by the FDA due to the concern of transference of prescription topicals to others. The Topi-Click directly addresses these concerns with needed innovation.

In addition, the audible clicker and tactile dosing feedback features of the Topi-Click allow patients to dispense a metered dose in a matter of seconds...even in complete darkness. The Topi-Click TDA protects the public and the drugs' success with features that help promote and attain compliance. With a metered dose comes more treatment success, thus also protecting profitability. Next, by combining patented dosing applicator technology with custom designed bottle shapes, pharmaceuticals are differentiated and pipelines are qualified. Custom Rx TDA devices to help because they help solve the growing topical dosing challenges with Topi-Click's laser focused innovation of features: 1. Applicator Part 2. Audible & Tactile Feedback of Metered Dosing; 3. Innovative Shape and Design for Recognition. See clips of highlights from Oprah, Dr. Phil and The Doctors that demonstrate its use with hormone prescriptions at www.Topi-Click.com

[Drug Deliveries Technology magazine reviews Topi-Click Innovation Award Nomination](#)

Drug Development & Delivery magazine provides information about the science, technology, and related business practices of drug delivery, drug development, pharmaceutical, biopharmaceutical, and specialty pharmaceuticals.

FIGURE 9
TECHNOLOGY INNOVATION NOMINEES

FIGURE 10
Individual Particle Electron Tomography for Structure Determination of a Highly Dynamic Protein

FIGURE 9 highlights the Topi-Click TDA as a Technology Innovation Nominee. It describes the device's features, including its metered dosing clicker, audible and tactile feedback, and its ability to protect the public and drug success. The article also mentions the device's use with hormone prescriptions and its recognition by Oprah, Dr. Phil, and The Doctors.

FIGURE 10 discusses the use of Individual Particle Electron Tomography for structure determination of a highly dynamic protein. It mentions the work of Dr. David Suck and his team at the University of Manchester, who used a cryo-electron microscope to determine the structure of a protein at the atomic level. The article also mentions the work of Dr. David Suck and his team at the University of Manchester, who used a cryo-electron microscope to determine the structure of a protein at the atomic level.

[Topi-Click use suggested in "Pharmacy Practice" textbook](#)

Fourth Edition: Mastering Community and Hospital Competencies

In [Pharmacy Practice](#), the chapter which discusses compounding and methods of dispensing features Topi-Click.

Pharmacy Practice for Technicians
Mastering Community and Hospital Competencies
Fourth Edition

Paradigm

Don A. Ballington
Robert J. Anderson

Nonsterile Pharmaceutical Compounding

8

Learning Objectives

- Define the proper techniques for weighing pharmaceutical ingredients, measuring, and transferring to a sterile container.
- Define the proper techniques for weighing pharmaceutical ingredients, measuring, and transferring to a sterile container.
- Define the proper techniques for weighing pharmaceutical ingredients, measuring, and transferring to a sterile container.

Selecting Medication Containers

With a much emphasis on chemical stability, it is important for the pharmacy technician to select the appropriate medication container to extend the beyond-use date, if such is possible. Standards for packaging the more common compounding are provided by USP in the Pharmacopeia (Pharmacopoeia).

Labeling and Cleanup

After the compounding operation, the preparation must be labeled with a precise quantity and the appropriate beyond-use date. The label should also include the name of the compound and the name of the pharmacist who prepared it. The label should also include the name of the compound and the name of the pharmacist who prepared it.

Visit www.Topi-Click.com to view full pages:



[Oprah Winfrey with Dr. Oz shows the Topi-Click being filled at a pharmacy and demonstrated by Dr. Christine Northrup and by Dr. Prudence Hall.](http://www.oprah.com/showinfo/The-Great-Debate-Should-You-Replace-Your-Hormones_1)

http://www.oprah.com/showinfo/The-Great-Debate-Should-You-Replace-Your-Hormones_1



- [Topi-Click helps raise almost \\$70,000 for the International Academy of Compounding Pharmacists Foundation at the 2011 IACP 6th Annual CCH Rx Mixer.](#)

Las Vegas came to the JW Marriott on Pennsylvania Avenue with a glittering array of celebrity impersonators to welcome the [2011 IACP 7th Annual Rx Mixer](#) — **Viva Las Vegas** benefitting the IACP Foundation. The casino-themed event began with cocktails at bars sponsored by [Topi-Click](#)..., and tempting dinner treats in the foyer that had guests standing in line. At the close of the evening, almost \$70,000 was raised to support the IACP Foundation’s research and educational programs.

- [Topi-Click helps raise over \\$90,000 for International Academy of Compounding Pharmacists Foundation at the 2010 IACP 6th Annual Compounders on Capital Hill Rx Mixer.](#)

The good times rolled for a good cause at the [2010 IACP 6th Annual Rx Mixer](#) – Mardi Gras. With over 277 guests in attendance, this year’s party was a hit! The event featured cocktails at bars sponsored by [Topi-Click](#)... More than \$90,000 in gross revenues from the auction was tallied-besting last year’s record.