

Maximize
the Benefits
of Topical
Medications



TOPI-CLICK®



SPECIFICATIONS

Each click (¼ turn of the Topi-Click® container base) moves the plunger (topical elevator) approximately 0.25ml of volume space, thereby dispensing 0.25ml of cream on average. The dispensing mechanism is similar to a syringe, but with more control of the plunger movement, so Topi-Click® is easier to use.

The Topi-Click® container can be used to dispense a total of approximately 5ml to 35ml of a topical preparation:

Topi-Click® is made of FDA grade plastic. UV block in plastic protects the contents from UV rays.

Topi-Click® can be labeled like a prescription vial.

Topi-Click® Benefits:

- Improves compliance for better outcomes of therapy through ease-of-use and convenience.
- Audible, visual and tactile clicks make it easy for even your most challenged patients to accurately dose medications.




- Demonstrations on Oprah, Dr. Phil and The Doctors national talk shows made Topi-Click® nationally recognized as the dosing applicator used for topicals.
- Distinguishable and elegant design that's unique from any other type of dispensing system to differentiate your treatments.
- Domed applicator pad allows medication to be applied without touching medication with the hands, thereby reducing the risk of spreading medicine to loved ones and pets.
- Easier to use than oral dosing syringes (doesn't resemble a needled syringe.)
- More pharmaceutically elegant than tubes or jars.
- Available in multiple colors (Blue, Pink and more) for differentiating various medications.
- Promotes additional compliance with a "Refill" reminder system so there are no gaps in treatments.
- User-friendly enough to be operated in the dark and almost impossible to discharge accidentally.
- Dispenses an average metered dose of 0.25 ml per click.

TOPI-CLICK®

makes it easy to accurately measure the dose and apply topical medications to the intended area, thereby decreasing the risk of improper dosing and contamination while improving compliance and outcomes of therapy.



Watch the video at www.Topi-Click.com to see the Topi-Click® in action and the highlights from Oprah, Dr. Phil, and The Doctors national talk shows.



Visit www.Topi-Click.com for more information